

	Date Issued: September 2017	Page No.: 1	Document No.: SOP-Biosafety-03
	Revision: 2.0	Subject: Biological Safety Cabinets	

1. Introduction and Scope

This SOP outlines the safety requirements and work practices to provide optimum contamination control and safety when working in a BSC.

When properly maintained and used in conjunction with good laboratory techniques, Biological Safety Cabinets (BSCs) provide effective primary containment for work with human and animal pathogens.

2. Scope

This procedure is applicable to all laboratories handling biohazardous materials classified as RG1/RG2 and RG2+ as defined in Queen's Biosafety Manual.

- In containment level 2 (CL2) laboratories, BSCs are used for procedures with the potential to produce infectious aerosols and for high concentrations or large volumes of infectious material.
- Most operations with RG2 material should be performed in a BSC unless a permit to work in an open bench has been approved by the Institutional Biohazard Committee (IBC)
- *All manipulations* of RG2 biohazardous agents that require level 2 containment with level 3 operational practices (commonly called CL2+) must be performed in a BSC.
- In containment level 1 (CL1) laboratories, BSCs can be used for some low risk procedures with RG1 materials. The procedures and verification requirements are recommended for the operation and maintenance of BSC.

3. Applicable Legislation, Standards, Guidelines

- Canadian Biosafety Standards, 3rd Edition, 2022 (CBS; Public Health Agency of Canada (PHAC) and Canadian Food Inspection Agency (CFIA))
- Canadian Biosafety Handbook, 1st Edition, 2016 (CBH; PHAC and CFIA)
- National Sanitation Foundation (NSF) Standard No. 49-2010 for the design, manufacture and testing of BSCs (NSF-49)

4. Responsibilities

4.1 Responsibilities of Department Heads

- Ensure that supervisors, employees and students are notified of their responsibilities for working with BSCs.

	Date Issued: September 2017	Page No.: 2	Document No.: SOP-Biosafety-03
	Revision: 2.0	Subject: Biological Safety Cabinets	

- Ensure that all employees and students have received instruction in the proper use of a BSC.
- Ensure that the components of this SOP and the applicable legislation are implemented in all facilities under the Head's authority.

4.2 Responsibilities of Supervisors

- Ensure that the correct type of BSC is used for the intended work.
- Contact EH&S to consult regarding purchase and placement in the laboratory.
- Ensure that BSC users receive safety training regarding best laboratory work practices and this SOP prior to commencing work with BSCs.
- Ensure that a BSC is decontaminated prior to being moved.
- Ensure that a BSC is tested and certified to NSF Standard 49 annually.
- Ensure that a BSC is tested and recertified after being moved. This is necessary because the HEPA filter can shift out of alignment or seals crack, and not function properly if the BSC is moved even a short distance.
- Ensure that no work is performed in a malfunctioning BSC.
- Receive report from EH&S regarding results of BSC testing.
- Coordinate maintenance, repair and new BSC equipment needs with EH&S.

4.3 Responsibilities of Users (Staff/Students)

- Understand how a BSC works.
- Ensure that the correct type of BSC is used for the intended work.
- Follow the procedures for BSC use described in this SOP.
- Notify the Laboratory Supervisor/Principal Investigator immediately in the event of BSC malfunction.
- Notify the Biosafety Officer if a BSC malfunction may have caused an exposure.

5. Definitions

5.1 Biological Safety Cabinet (BSC)

Biosafety cabinets (BSCs) are essential engineering devices for protecting the safety of lab personnel, the environment, and samples when working with biohazardous agents. BSCs use a combination of directional airflow and a High Efficiency Particulate Air (HEPA) to filter and remove bacteria, spores and viruses with 99% efficiency.

The unique directional airflow in a BSC consists of:

- (1) HEPA-filtered air that passes over the work surface to provide a clean working environment;
- (2) an "air curtain" at the main opening that prevents the mixing of potentially

	Date Issued: September 2017	Page No.: 3	Document No.: SOP-Biosafety-03
	Revision: 2.0	Subject: Biological Safety Cabinets	

- contaminated air from the laboratory with the HEPA-filtered air in the BSC;
- (3) continuous negative pressure in the BSC that is created by a continual draw of air into the BSC. This prevents the air inside from escaping into the laboratory;
 - (4) the air that is exhausted out of the cabinet into the laboratory is HEPA-filtered to prevent the release of particulate contaminants.

The type of BSC determines whether it provides only protection for the operator, or protection for both the operator and the material in use i.e. whether or not it provides an aseptic environment.

Note that there are also “laminar flow cabinets/hoods/clean air hoods” that provide only product protection, but no operator protection. These are not BSCs (although they may look somewhat similar), so it is important to know the difference and be sure that you are using a BSC when operator protection is needed. Most often these have air blowing towards the operator, which is not the case with a BSC.

There are three classes of BSC and selection of the proper class requires careful evaluation of the activities to be carried out. A table and figures detailing the features of the most common BSCs at Queen’s are provided in the Appendix I.

Class I - have un-recirculated airflow away from the operator that is discharged to the atmosphere after filtration through a HEPA filter. They provide good operator protection but no product protection.

Class II types A1, A2, B1, and B2 - are designed for personnel, product and environmental protection. The main difference between the types are the ratio of air exhausted from the BSC to that which is recirculated within the BSC, and the type of exhaust system present. Type A cabinets usually exhaust filtered cabinet air into the room, but are sometimes ducted out of the building by means of a "thimble connection".

Type A1 cabinets are never hard-ducted (i.e. never exhausted by a hard duct to the outside of the building) and are not suitable for work with volatile toxic chemicals or volatile radionuclides because 70% of the air is recirculated which could cause a dangerous buildup of the toxic materials inside the BSC, or inside the containment zone. **These cabinets require a special test that is no longer performed by Con-Test. This means that these types of cabinets can no longer be used for PHAC regulated biohazardous materials. PHAC regulated material includes risk group (RG)2 human & animal pathogens and toxins. It also includes viral vectors and cell lines listed as RG2 on the ePathogen/registry database.**

	Date Issued: September 2017	Page No.: 4	Document No.: SOP-Biosafety-03
	Revision: 2.0	Subject: Biological Safety Cabinets	

Type A2 cabinets, appear similar to type A1 cabinets; however, they have a greater inflow velocity and always have negatively pressured contaminated plenums or positively pressured contaminated ducts/plenums surrounded by negatively pressured ducts/plenums. This design feature ensures that the potential leaks in the positively pressured ducts or plenums are drawn inward rather than out into the containment zone. This type of BSC is suitable for work with minute amounts of volatile toxic chemicals and radionuclides, if ducted out of the building with a thimble connection.

Type B cabinets are hard-ducted through a dedicated duct exhausted to the atmosphere after passage through a HEPA filter.

- Type B1 is suitable for work with low levels of volatile toxic chemicals and trace amounts of radionuclides, as long as the work is performed towards the back of the cabinet where air is discharged directly to the outside atmosphere.
- Type B2 is suitable for work with volatile toxic chemicals and radionuclides because they are hard-ducted and the airflow within them is designed so that air is never recirculated within the BSC or within the containment zone.
- Reversal of airflow from the face of a BSC, also known as a **puff-back**, can occur in Class II type B2 BSCs, for example upon failure of the HVAC system, power, or the exhaust fan serving the BSC. Every effort is to be made to address puff-backs mechanically (CBS Matrix 3.7).
- Consideration should also be given to the amount of air required to operate this type of cabinet as it may lead to additional adjustments to balance the airflow in the containment zone.

Type C1 cabinets is a new type of BSC that can be used in a recirculating Type A-mode for standard microbiological work, or can be connected to an exhaust system to function in Type-B mode for handling hazardous chemical vapors or radionuclides. There are several advantages to this design over a traditional Type B BSC, so they should be considered if direct exhaust is required.

Class III – are totally enclosed and gas-tight with HEPA filtered supply and exhaust air. Work is performed with attached long-sleeved gloves. They are designed for work with level 4 pathogens. None are in use at Queen's University.

Clean Air Benches

Clean air benches do NOT protect the worker and are NOT BSCs. Clean air benches have HEPA filtered laminar airflow *towards the worker*. They provide a flow of clean air over the product and protect it from contaminants in the environment. A clean air bench provides

	Date Issued: September 2017	Page No.: 5	Document No.: SOP-Biosafety-03
	Revision: 2.0	Subject: Biological Safety Cabinets	

product protection only. The worker is directly exposed to aerosols and particulates from the work. Clean air benches are not to be used for work with biohazard risk group 2 material, hazardous volatile chemicals (or particulates), or radioisotopes.

HEPA filter

HEPA filters can remove at least 99.97% of airborne particles 0.3 μm in diameter. Particles of this size are the most difficult to filter and are thus considered the *most penetrating particle size* (MPPS). *Particles that are larger or smaller than 0.3 μm are filtered with even higher efficiency.* However, **HEPA filters do not prevent gases from passing through.**

HEPA filters are composed of a mat of randomly arranged fibres. Particles are trapped (stick to a fibre) by one of the following three mechanisms:

- 1) *Interception*, where particles following a line of flow in the air stream come within one radius of a fibre and adhere to it.
- 2) *Impaction*, where larger particles are unable to avoid fibres by following the curving contours of the air stream and are forced to embed in one of them directly; this effect increases with diminishing fibre separation and higher air flow velocity.
- 3) *Diffusion*, an enhancing mechanism is a result of the collision with gas molecules by the smallest particles, especially those below 0.1 μm in diameter, which are thereby impeded and delayed in their path through the filter; this behaviour raises the probability that a particle will be stopped by either of the two mechanisms above; it becomes dominant at lower air flow velocities. Diffusion predominates below the 0.1 μm diameter particle size. Impaction and interception predominate above 0.4 μm . In between, near the 0.3 μm MPPS, diffusion and interception predominate.

Mechanical Operation of a Biological Safety Cabinet

- **Operation Panel and Main Features:** The functions of a BSC are controlled through a panel located on the front of the cabinet above or beside the sash (window). Switches on the panel control the lighting, the blower/fan system, the electrical outlets, and the alarm system bypass.
- **Magnehelic Gauge:** Routinely monitor the pressure on the magnehelic gauge so that you will know what it should be for your cabinet since normal reading differs from one cabinet to another. The magnehelic gauge provides a "gross" or approximate indication of HEPA filter loading, i.e., the volume of particulate matter the filter has accumulated as the cabinet operates. It provides one indication of how rapidly the filter's capacity is being diminished. This is measured by reading either motor blower suction or pressure.

	Date Issued: September 2017	Page No.: 6	Document No.: SOP-Biosafety-03
	Revision: 2.0	Subject: Biological Safety Cabinets	

- Blower:** The blower motor should be turned on at least 5 minutes prior to using the BSC. To maintain the proper airflow pattern, the sash (window) should only be opened to the height marked on the front panel (usually a maximum of 10 inches). Although the BSC may be turned off after decontamination as described under "Upon completion of work" (page 6), many users choose to leave the BSC running continuously to maintain a clean environment at all times. In addition, if a BSC is used many times during a day, it should be left running until all work is complete because frequent turning on and off of the BSC decreases the life expectancy of the motor.
- Alarm:** There should be an alarm that will sound if the air circulation in the cabinet has been compromised (usually because the sash or window has been raised too far, but also for other reasons such as a blocked HEPA filter). A BSC must not be used when an alarm sounds or a warning light is on. An alarm must never be disengaged as it indicates improper airflow which affects the cabinet's performance and may endanger the researcher or the experiment.
- On older BSCs that do not have alarms, it is particularly important to routinely check the magnehelic gauge to ensure that the air flow pressure remains constant.
- Lighting:** There are two types of lighting inside the BSC: fluorescent and Ultraviolet (UV). The fluorescent light should normally be ON while work is being done within the cabinet. Although the direct exposure to UV light (approximately 260 nm wavelength) can reduce the number of pathogenic microorganisms on exposed surfaces and in air, **the routine use of UV lamps to decontaminate a BSC is not recommended for the following reasons:**
 - UV light has poor penetrating power.
 - UV irradiation is ineffective if a
 - microorganisms is protected by dust, dirt, or organic matter. A liquid chemical disinfectant should be the primary method of cleaning and disinfecting the interior of a BSC.
 - UV irradiation does not penetrate into cracks or through the grilles of a BSC.
 - UV irradiation can cause deterioration of various materials, including certain plastics and tubing
 - the accumulation of dust, dirt, grease or clumps of microorganisms reduce its
 - germicidal effects;
 - exposure to UV light is hazardous: it may result in severe eye damage and burns to the skin.

Nevertheless, UV lights may be **useful in certain situations, if properly maintained,**

	Date Issued: September 2017	Page No.: 7	Document No.: SOP-Biosafety-03
	Revision: 2.0	Subject: Biological Safety Cabinets	

such as when using spore-forming bacteria since bacterial spores are resistant to chemical disinfection.

- UV irradiation of the work area should only be used as a secondary method of disinfection in the cabinet. Never rely on UV irradiation alone to disinfect a contaminated work area.
 - UV lamps **must** be turned off whenever the laboratory is occupied unless the sash can be lowered so that it is completely closed.
 - The room should never be occupied when UV lamp is in use (UV light will damage the human eye and skin very quickly).
 - The vertical sliding window must be in the closed position. If it has a hinged window, use manufactures enclosure panel.
 - The lamp must be wiped free of dust/linT weekly.
 - The UV wavelength must be checked periodically.
- **Outlets:** Electrical devices can be used inside the BSC by using the plugs provided.

6. Safety precautions

- Only one person should work in a BSC at one time (even in 6 foot long BSCs). Two people will disrupt the airflow, creating turbulence that could result in biohazardous aerosols coming out of the cabinet, or cause cross-contamination of work in the cabinet. i.e. it can be hazardous for you and/or your science.
- Large BSCs are appropriate for experiments that need lots of space, not for two people to work at one time.
- Wear protective gloves that cover the cuffs of lab coat sleeves to prevent contaminated air from entering the sleeve. Lab coats (or closed front gowns) with fitted cuffs rather than loose sleeves are recommended.
- Equipment creating air movement (e.g., vacuum pumps, centrifuges) may affect the integrity of the airflow and should not be used within the BSC.
- To help protect the fragile air curtain at the front of the cabinet, BSCs must be located as far as possible away from areas where airflow patterns may be disrupted e.g., room air supply and exhaust grilles, doors, open windows, high traffic areas, and large pieces of equipment that generate heat (CBS Matrix 3.7.6).
 - Consult EH&S if your BSC appears to be located in an area where the air curtain may be disrupted. EH&S can do a smoke test on the air curtain, and propose

	Date Issued: September 2017	Page No.: 8	Document No.: SOP-Biosafety-03
	Revision: 2.0	Subject: Biological Safety Cabinets	

remedial actions.

- Windows that open should be kept closed when the BSC is in use.
- **An open flame should not be used in a BSC.** Natural gas and propane should not be used in a BSC; sustained open flames (e.g., Bunsen burner) in BSCs are prohibited. On-demand open flames (e.g., touch-plate microburners) are to be avoided as they create turbulence in the BSC, disrupt airflow patterns, and can damage the HEPA filter (CBS Matrix 4.6). **Non-flame alternatives (e.g., microincinerator, or sterile disposable inoculation loops) should be used whenever possible.**
- The HEPA filters in the BSCs remove particulates from air, but they are **not effective at collecting chemical gases or vapours**. If you need to use such material in a BSC, contact the Biosafety Officer for advice.
- **Routine use of UV lamps to decontaminate a BSC is not recommended**
- **DO NOT use the BSC if the ALARM sounds or if there are other indications of cabinet malfunction such as no airflow, reduced pressure on manometric gauge (drop > 0.2), or unusual noises.**

7. Work practices

The work practices indicated in this procedure are for general recommendations. Always ensure you are also familiar with the cabinet specifications and manufacturer operational and verification requirements.

7.1 BSC Start up

1. Turn **OFF** the UV lamp (if it is on) and turn **ON** the fluorescent light.
2. Turn **ON** the blower (if it is off; many labs leave the blower running all the time) and allow it to run for at least 5 minutes before starting work in the cabinet.
3. Let the cabinet run for at least 5 minutes to purge its content and equilibrate.
4. **Check BSC flow is within operational range.** If the differential pressure drops by more than 0.2 then do not use the cabinet and call for service.
5. **Check that the sash is at the appropriate height.** Adjust stool height so that the user's underarms are level with the bottom of the sash.
6. If present, test the **airflow alarm is functioning** and ensure it is switched to the "on" position.
7. Fill out the BSC maintenance log sheet by initializing all the applicable steps .

7.2 BSC setup

1. Disinfect the interior surfaces with a disinfectant appropriate for the biohazardous agents

	Date Issued: September 2017	Page No.: 9	Document No.: SOP-Biosafety-03
	Revision: 2.0	Subject: Biological Safety Cabinets	

used in the laboratory (e.g. 70% Ethanol (EtOH), Virox, an H₂O₂ disinfectant). Note that bleach can corrode the stainless steel and make the BSC difficult to decontaminate so it should not be used for routine decontamination. If bleach or other corrosive disinfectant is required (e.g. for spill clean-up) then be sure to thoroughly rinse the surface.

2. Assemble all materials required for manipulation and load into the BSC.
 - a. Care should be taken not to overcrowd or block the front or rear grilles to prevent the appropriate airflow patterns from being compromised.
 - b. Whenever possible, place all materials needed for a procedure inside the cabinet before starting.
 - c. Place supplies, equipment and papers well back from the front of the cabinet, positioned so that air intake or exhaust grills are not obstructed.
 - d. Never put anything on the grill at the front opening of the cabinet.
 - e. Do not block the air openings/grill at the back of the cabinet.
 - f. Segregate non-contaminated ("clean") items from contaminated ("dirty") items. Work should always flow from "clean" to "dirty" areas
 - g. **Disinfectant effective against the agents in use should be placed inside the BSC** so that you do not have to reach out with contaminated gloves to get it.
 - h. Place **aerosol generating equipment** (e.g., vortex mixer, sonicator) towards the back of the BSC, without blocking the rear grille.
 - i. When there is significant **potential for splatter or splashes** to occur during manipulations of infectious material or toxins, the work area should be lined with a plastic-backed absorbent pad.



3. After loading material in the BSC, allow sufficient time for the air to purge and the airflow to stabilize before initiating work. This will be specified in the manufacturer's instructions, and is generally 3-5 minutes.

	Date Issued: September 2017	Page No.: 10	Document No.: SOP-Biosafety-03
	Revision: 2.0	Subject: Biological Safety Cabinets	

7.3 Working in the BSC

1. Perform operations as far to the rear of the work area as reasonable. Ensure that elbows and arms do not rest on the grille or work surface.
2. Move arms slowly and move straight out of the cabinet perpendicular to the front opening; do not sweep arms across the front of the cabinet. Do not walk quickly in front of a cabinet when someone else is working. Movement of arms into and out of the cabinet can disrupt airflow, which can allow contaminants to enter or escape the BSC.
3. Routinely monitor the pressure on the pressure (magnehelic) gauge so that you will know what it should be for your cabinet since the normal reading differs from one cabinet to another.
4. Material should be discarded in a **waste container located towards the rear of the cabinet workspace. Do not discard materials in containers outside of the cabinet.**
5. If using a vacuum aspirator, please refer to SOP-Biosafety-01 for proper set up use and maintenance.

7.4 Completion of work:

1. Upon completion of work, allow sufficient time for the air in the BSC to purge (i.e., pass through the filter) before disrupting the air curtain by removing hands or unloading material from the BSC. The purge time will vary by model and can be up to several minutes.
2. Surface-disinfect objects before removal from the cabinet. Remember that aerosols generated during operations in the cabinet such as pipetting might have contaminated objects in the cabinet, so there does not have to have been a spill for this step to be necessary.
3. Remove all containers and equipment from the cabinet
4. Disinfect the interior surfaces of the BSC, including sides, back, lights, and interior of the glass, with a disinfectant effective against the pathogens in use, allowing an appropriate contact time. If a corrosive disinfectant is used, the

	Date Issued: September 2017	Page No.: 11	Document No.: SOP-Biosafety-03
	Revision: 2.0	Subject: Biological Safety Cabinets	

surface should be thoroughly rinsed with water after disinfection to avoid corrosion of the stainless steel surfaces. If you are using a chloride type disinfectant, wipe down interior surfaces with a 70% alcohol solution. This will protect stainless steel interiors from corrosion.

5. If necessary, the cabinet should also be monitored for radioactivity and decontaminated when necessary.
6. If appropriate turn off the blower and fluorescent lamp, and turn on the UV light if appropriate. UV lamps are neither recommended, nor required in Biological Safety Cabinets; however, if you are going to use one please follow the safety guides below to ensure you and your lab's safety.
7. Remove gloves and gowns and wash your hands as the final step in safe microbiological practices

7.5 Spills inside the cabinet

1. Clean up spills as soon as they occur. Follow the steps outlined in your Laboratory Emergency Procedures – Biological Material Spills.
2. Notify other lab members of the spill and notify your supervisor after cleaning the spill.
3. Submit a spill report using the [Spill Report](#) form from the EH&S website.

7.6 Instrument failure

If alarm or other indication of failure happens while using the cabinet:

1. Remain Calm.
2. Immediately stop work.
3. Close all primary containers.
 - If not working with biohazardous materials, surface decontaminate items, remove the, from the BSC, post "Do Not Use" signage and report problem to supervisor; if working with biohazardous materials continue to follow the instructions below.
4. Remove gloves and discard in BSC waste.
5. Turn off BSC and close sash.
6. Wash hands.
7. Post "Do Not Use" signage.
8. After 30 minutes, don clean PPE, surface decontaminate and remove any biohazardous material.

	Date Issued: September 2017	Page No.: 12	Document No.: SOP-Biosafety-03
	Revision: 2.0	Subject: Biological Safety Cabinets	

9. Decontaminate work surface, close sash, leave "Do Not Use" signage on BSC..
10. Report problem to Supervisor.

7.8 Maintenance/Certification of Biological Safety Cabinets requirements

Maintenance and verification requirements depend on the classification of the cabinet. Use the table below to know the requirements for your BSC.

BSC type	Class I	Class II A1	Class II A2	Class II B1	Class II B2
NSF/ANSI 49 periodic certification	Recommended annually	Annually*	Annually*	Annually*	Annually*
Puff back test verification	Not required	Not required	If thimble connected After Installation modification, system failure, or repair	Not required	After Installation modification, system failure, or repair
Regular Maintenance	Recommended monthly	6 months	6 months	6 months	6 months
Decontamination	Before relocation or instrument/laboratory decommissioning				

***NSF certified biosafety cabinetry should be used for level 2 containment.**

- NSF/ANSI Standard 49, includes basic requirements for design, construction and performance to provide personnel, product and environmental protection; reliable operation; durability; cleanability; noise level and illumination control; vibration control; and electrical safety. In addition, the standard includes detailed test procedures and informational annexes, including recommendations for installation, field certification tests, and decontamination procedures.
- When used for biosafety level 2 containment, biological safety cabinet performance must be tested and certified (as specified in CBS Matrix 5.1):
 - upon initial installation in the laboratory
 - annually thereafter
 - when moved from one area to another within the same room, or from one room to another

	Date Issued: September 2017	Page No.: 13	Document No.: SOP-Biosafety-03
	Revision: 2.0	Subject: Biological Safety Cabinets	

- whenever maintenance is carried out on internal parts, and whenever filters are changed
- Where a custom enclosure or the design of a BSC or similar enclosure (e.g. Biobubble) does not permit certification in accordance with NSF/ANSI 49, it is to be verified to meet the manufacturer's specifications, with minimum parameter values specified in Matrix 5.1 of the CBS

7.9 Regular maintenance procedures

1. BSC regular maintenance must be performed monthly.
2. Use the link [BSC MAINTENANCE LOG.pdf](#) located on the EHS website to access the form and it is also found at the end of this SOP.

7.9.1 Cleaning Work Surface

1. Put on lab coat, gloves and safety glasses.
2. Make sure BSC is empty of any items, equipment, or waste containers.
3. Allow the cabinet to operate for 5 minutes with no activity, which should purge airborne contaminants from the work area.
4. Wet wipe the top surface of the outside edge of the BSC and front grille a total of 3 times with the pre-soaked paper towels with 1% sodium hypochlorite.
5. Place used paper towels into a regular waste bag.
6. Decontaminate work surface before removing it from the cabinet by spraying/wiping with 70% ethanol.
7. Remove work surface: **consult the operation manual of your cabinet and ensure this surface is removed safely.**
8. Thoroughly spray the underside of the work surface with 70% ethanol. Let sit for 5 minutes.
9. Assess the underside of the work surface and clean any gross contamination present by wet-wiping 3x the entire surface with paper towels pre-soaked with decontamination solution.
10. Rinse by wet-wiping with paper towels pre-soaked with tap water if necessary. Using tongs or another long reaching tool is recommended to clean hard to reach surfaces.
11. Collect all paper towels into the regular waste bag.

7.9.2 Cleaning the Front Grille

1. Remove the front grille: **consult the operation manual of your cabinet and ensure this part is removed safely.**

	Date Issued: September 2017	Page No.: 14	Document No.: SOP-Biosafety-03
	Revision: 2.0	Subject: Biological Safety Cabinets	

2. Wet-wipe the underside of the grille for a total of 3 times with pre-soaked paper towels.
3. Place used paper towels into a regular waste bag.
4. Scrub any gross contamination remaining on the grille with appropriate scrubbing tool. Remove all loosened debris by wet-wiping another three times.
5. Spray/rinse with water to remove sodium hypochlorite residue.
6. Spray top of grille and underside of grille thoroughly with 70% ethanol. Let sit for 5 minutes.

7.9.3 Cleaning BSC Catchbasin

1. Using a flashlight and an extendable mirror, examine the extent of cleaning required. Note any sharps or broken glass.
2. Saturate the entire surface of the catch basin by misting with decontamination solution to prevent fly away debris. You want to capture as much solid material and not let it get sucked up into the HEPA filters.
3. Let sit for 5 minutes, keep spraying, do not let it dry out.
4. Using an appropriate tool (such as: tongs, plastic scraper, mini dustpan, toilet brush or long handled scrub brush) remove as much loose solid material, sharps or broken glass. Place sharps or broken glass into a sharps container.
5. Collect debris into a plastic bag inside the BSC if you are positive no sharps are present. If the debris contains sharps that you cannot pick out, place the wet, sharps-containing muck into a wide mouth sharps container.
6. For dried on gross contamination, scrub catchbasin with a plastic scrub brush. A toilet brush is ideal because of the long handle.
7. Place paper towels soaked in decontamination solution in the catchbasin. Make sure to cover entire surface.
8. Let sit for 15 minutes or the recommended contact time for the decontamination solution. Remove paper towels and dispose of in regular waste bags.
9. Repeat steps 6 to 8 to ensure all debris is removed.
10. Rinse catchbasin with tap water three times. Failure to thoroughly rinse corrosive chemicals will result in damage to the entire catchbasin.
11. Spray catchbasin with 70% ethanol and let it air dry.

7.9.4 Cleaning the Rest of the BSC

1. Replace front grille

	Date Issued: September 2017	Page No.: 15	Document No.: SOP-Biosafety-03
	Revision: 2.0	Subject: Biological Safety Cabinets	

2. Replace work surface
3. Surface disinfect using 70% ethanol to decontaminate all interior work surfaces, making sure to clean the back and side of the cabinet (sidewalls, back wall), inside of sash (interior of the glass), and work surface.
4. If using a chloride type disinfectant, after contact time, wipe down interior surfaces with a 70% alcohol solution to protect stainless steel interiors from corrosion.
5. The drain pan should be emptied into a collection vessel containing disinfectant; and the drain valve can be disinfected using a flexible tube.
6. Any spilled fluid and disinfectant solution on the work surface should be absorbed with paper towels and discarded into a biohazard bag.
7. Hands should be washed whenever gloves are changed or removed
8. If necessary, the cabinet should also be monitored for radioactivity and decontaminated when necessary.

7.9 Cont.. After completing the tasks

3. Remove and dispose of gloves.
4. Remove and disinfect the lab coat and safety glasses using appropriate disinfectants and procedures.
5. Record the date in the BSC maintenance log by adding your initials to the corresponding date of the month.

7.10 Information and Enquires

To arrange BSC maintenance or testing and certification please email tom.martinek@queensu.ca. For technical questions or maintenance issues contact the BSC technician (ext. 74976); for other questions contact the University Biosafety Officer (ext. 77077).

8. Training requirements

Users of this procedure must be trained on

- Biosafety Level 1 and 2
- WHMIS 2015
- WHMIS annual refresher
- EHS- Fundamentals of working safely in a Class II Biological Safety Cabinet

9. References

[Biological Safety Cabinet Decontamination: 10 Step Guide \[Video\]](#)

[biosafety-cabinet-preventative-maintenance.pdf](#)

[66340 LM Biosafety ebook AV V6.pdf](#)

[Best Practices for Biosafety Cabinet Maintenance and Use | Lab Manager](#)

	Date Issued: September 2017	Page No.: 16	Document No.: SOP-Biosafety-03
	Revision: 2.0	Subject: Biological Safety Cabinets	

[LM-SOP-BSC Routine Maintenance.pdf](#)

10. **Approval and Revision History:** approved by the University Biohazards Committee

Revision History

Date	Revised by	Changes
V1.2009	Shelagh Mirski	Initial release
V2. 2017	Shelagh Mirski	Update to include information from CBS 2nd Edition 2015 and CBH 1st Edition, 2016
V3. 2023	Raico Laria Lamela & Natalie Roy	Updating what
V4. 2026	Raico Laria Lamela & Natalie Roy	Update CBS info, add revision history & add maintenance log information. Maintenance procedure added. Document reformatted to EHS template SOP.

	Date Issued: September 2017	Page No.: 17	Document No.: SOP-Biosafety-03
	Revision: 2.0	Subject: Biological Safety Cabinets	

Appendix I - Queen's University SOP-Biosafety-03 Biological Safety Cabinets

Table and Figures showing different types of BSC, copied from the Canadian Biosafety Handbook, 1st Edition, 2016, Chapter 11


Table 11-1: Summary Table of Key Characteristics of Class II Biological Safety Cabinets (BSCs)

	Type A1	Type A2	Type B1	Type B2
Minimum average inflow velocity through front opening	0.38 m/s [75 fpm]	0.51 m/s [100 fpm]	0.51 m/s [100 fpm]	0.51 m/s [100 fpm]
Air patterns	30% of the air is exhausted out of the BSC and 70% of the air is recirculated within the BSC	30% of the air is exhausted out of the BSC and 70% of the air is recirculated within the BSC	>50% of the air is exhausted out of the BSC and <50% of the air is recirculated within the BSC	100% of the air is exhausted out of the BSC
HEPA-filtered downflow air	Composed of mixed downflow and inflow from common plenum	Composed of mixed downflow and inflow from common plenum	Inflow air	Drawn from the containment zone or from the outside atmosphere
HEPA-filtered exhaust air	Recirculated to the containment zone or directly to the outside atmosphere	Recirculated to the containment zone or directly to the outside atmosphere	Exhausted through dedicated exhaust plenum to the outside atmosphere	Exhausted through dedicated exhaust plenum to the outside atmosphere
Type of exhaust	Can be thimble connected	Can be thimble connected	Hard-ducted	Hard-ducted
Contaminated ducts and plenums	Negatively pressured or surrounded by negatively pressured ducts or plenums; plenum may be positively pressured in some models	Negatively pressured or surrounded by negatively pressured ducts or plenums	Negatively pressured or surrounded by negatively pressured ducts or plenums	Negatively pressured or surrounded by negatively pressured ducts or plenums
Work with volatile toxic chemicals and radionuclides	No	Minute amounts if exhausted through thimble connection if present	Low levels of volatile toxic chemicals and trace amounts of radionuclides	Yes

	Date Issued: September 2017	Page No.: 18	Document No.: SOP-Biosafety-03
	Revision: 2.0	Subject: Biological Safety Cabinets	

Figure 11-2: Illustration of a Class II Type A1 Biological Safety Cabinet (BSC) (with a Positively Pressured Contaminated Plenum)


Cabinet exhaust may be recirculated into the room or vented to the outside atmosphere through an air gap type (thimble) connection, as shown. Purple shading indicates positively pressured contaminated plenum.



	Date Issued: September 2017	Page No.: 19	Document No.: SOP-Biosafety-03
	Revision: 2.0	Subject: Biological Safety Cabinets	

Figure 11-3: Illustration of a Class II Type A2 Biological Safety Cabinet (BSC)


Cabinet exhaust may be recirculated into the room or vented to the outside atmosphere through an air gap type (thimble) connection, as shown. Cabinet shown has a negatively pressured plenum.



	Date Issued: September 2017	Page No.: 20	Document No.: SOP-Biosafety-03
	Revision: 2.0	Subject: Biological Safety Cabinets	

Figure 11-5: Illustration of a Class II Type B2 Biological Safety Cabinet (BSC)

Cabinet is vented to the outside atmosphere through a hard-ducted connection, as shown.





Biological Safety Cabinet
MAINTENANCE LOG

Laboratory Building & Room Number _____ Month _____ Year _____

BSC Model # _____ BSC Serial Number# _____

DURING USE	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31
Check BSC in operational range (if applicable)																															
Check sash height																															
Before work begins, disinfect cabinet and sash																															
After work is completed, disinfect cabinet & sash and surface disinfect objects																															
Check drain valve under cabinet is closed																															
Check intake air grille not blocked																															
Let cabinet run for 5 mins after work is complete, to purge atmosphere before turned off.																															
MONTHLY																															
Remove & disinfect grille area																															
Clean catch basin																															

Person responsible for BSC maintenance: _____