



Aspiration MAINTENANCE LOG

Laboratory Building & Room Number _____

Month _____

Year _____

**Aspiration in-line filters typically are replaced once resistance occurs, or mold growth occurs.*

BEFORE USE	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31
Visual inspection of tubing and in-line filter																															
Check aspiration flask liquid level																															
WEEKLY																															
Replace bleach in aspiration flask (after 7 days)																															
AS REQUIRED																															
Replace in-line filter																															

Person responsible for aspiration maintenance: _____