

CENTRIFUGE MAINTENANCE LOG

Laboratory Building & Room Number _____ Month _____ Year _____

Centrifuge Model # _____ Centrifuge Serial Number# _____

**Use Manufacturer's recommendations for your centrifuge model when performing maintenance check.*

***Initial maintenance items only when applicable*

BEFORE USE	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31
Visual inspection of rotors, canisters, adapters and caps																															
Lid hydraulics, operation check																															
MONTHLY																															
Clean inside centrifuge																															
Clean outside centrifuge																															
Clean rotor, adapters, caps, o-rings																															
Coat o-rings with vacuum grease																															
Grease bucket pivot/pins/sockets																															
Wash filters																															
Change balance tubes																															

Person responsible for centrifuge maintenance: _____